

Electrochemical Testing and consolation Charges-2026

(@All prices Excluding 18% GST)

| Sl. No | Technique | Rs | Description/ Specifications |
|--------|--|-------|-----------------------------------|
| 1 | Linear Sweep Voltametry (LSV) | ₹2999 | per electrode |
| 2 | Cyclic Voltammety (CV) | ₹2999 | per electrode |
| 3 | Electrochemical Impedance Spectroscopy (EIS) | ₹2999 | per electrode |
| 4 | Chronopotentiometry (CP) OR Galvonostatic Charge-Discharge (GCD) | ₹2999 | per electrode |
| 6 | Differential Pulse Voltammety (DPV) | ₹2999 | per electrode |
| 7 | Galvonostatic Intermittent Titration Technique (GITT) | ₹2999 | per electrode |
| 8 | HER Studies (LSV, CV, EIS) | ₹7999 | per sample |
| 9 | OER Studies (LSV, CV, EIS) | ₹7999 | Per sample |
| 10 | Stability Studies | ₹2999 | Up to 12h after that Extra ₹299/h |
| 11 | Electrode Preparation/Utilization, (Ni foam, Torey carbon, Pt/C, Reference Electrode, Electrolyte-KOH, H ₂ SO ₄ , Nafion, Binder, etc) | ₹1999 | per electrode, per sample |
| 12 | Supercapacitor (CV, GCD, EIS)-3 electrode (Half Cell) | ₹9999 | per electrode |
| 13 | Supercapacitor (CV, GCD, EIS)-2 electrode (Full cell) | ₹7999 | per Cell |
| 14 | Coin Cell Fabrication | ₹699 | per Cell |
| 15 | Supercapacitor (CV, GCD, EIS)-Coin Cell | ₹5999 | per Cell |

Electrochemical Sensor (Biomolecules-Glucose, Uric Acid, Ascorbic Acid, Creatinine)

| | | | |
|----|--|-------|---------------|
| 16 | Cyclic Voltammety (CV) | ₹2999 | per electrode |
| 17 | Electrochemical Impedance Spectroscopy (EIS) | ₹2999 | per electrode |
| 18 | Differential Pulse Voltammety (DPV) | ₹2999 | per electrode |
| 19 | Interference, stability, repeatability | ₹4999 | per electrode |

Electrochemical Results Analysis & Reporting Charges

| Sl. No | Technique | Rs | |
|--------|-----------|----|--|
|--------|-----------|----|--|

| 1 | HER Studies (LSV, CV, EIS, Tafel, stability) | ₹14999 | ** |
|--|--|--------|----|
| <p>Analysis & Reporting Includes: **</p> <p>LSV Data Processing & Interpretation</p> <ul style="list-style-type: none"> Plotting raw and corrected LSV curves with proper baseline adjustment Determining onset potential and overpotential at benchmark current densities (e.g., η_{10}) Identifying activation-controlled and diffusion-controlled regions <p>CV Data Analysis</p> <ul style="list-style-type: none"> Plotting CV curves with scan-rate corrections for reliable comparison Identifying redox features and catalytically active sites Extracting double-layer capacitance (C_{dl}) and estimating ECSA for intrinsic activity comparison <p>EIS Data Modelling</p> <ul style="list-style-type: none"> Generating clear Nyquist and Bode plots Fitting data using appropriate equivalent circuit models Extracting R_s, R_{ct}, C_{dl}, and diffusion parameters to understand charge-transfer kinetics Interpreting interfacial behaviour and catalyst–electrolyte interactions <p>Tafel Plotting & Mechanistic Insights</p> <ul style="list-style-type: none"> Plotting η vs $\log j$ with proper linear fitting Calculating Tafel slope and exchange current density (j_0) Identifying the dominant HER reaction pathway (Volmer, Heyrovsky, or Tafel mechanism) <p>Stability Study Evaluation</p> <ul style="list-style-type: none"> Plotting CA/CP curves and analysing performance retention over time Comparing pre- and post-stability LSV to assess durability Identifying degradation trends and possible structural/chemical changes <p>Comprehensive Reporting</p> <ul style="list-style-type: none"> Providing refined, publication-ready plots (LSV, CV, EIS, Tafel, CA/CP) Presenting extracted parameters in structured tables Offering clear interpretation of catalytic behaviour and performance trends Delivering a concise conclusion summarising key outcomes <p>Importance of These Analyses for HER Research & Manuscripts</p> <ul style="list-style-type: none"> These analyses collectively establish activity, kinetics, mechanism, and stability, which are the four mandatory pillars of any high-quality HER study. Together, they provide all core electrochemical evidence required for writing a complete, publication-ready manuscript, enabling comparison with state-of-the-art catalysts and satisfying reviewer expectations for depth, reproducibility, and scientific rigour. | | | |

| | | | |
|---|---|--------|----|
| | | | |
| 2 | OER Studies (LSV, CV, EIS, Tafel, stability) | ₹14999 | ** |
| | <p>Analysis & Reporting Includes: **</p> <p>LSV Data Processing & Interpretation</p> <ul style="list-style-type: none"> Plotting raw and corrected LSV curves with baseline adjustment Determining onset potential and overpotential at standard current densities (e.g., η_{10} / η_{20} / η_{50}) Identifying kinetic vs. diffusion-limited regions Smoothing and noise reduction for high-quality comparative plots <p>CV Data Analysis</p> <ul style="list-style-type: none"> Plotting CV curves at multiple scan rates for reliable electrochemical profiling Identifying redox transitions associated with OER-active species (e.g., $\text{Ni}^{2+}/\text{Ni}^{3+}$, $\text{Co}^{2+}/\text{Co}^{3+}$) Calculating double-layer capacitance (C_{dl}) and estimating ECSA for normalization of intrinsic activity <p>EIS Data Modelling</p> <ul style="list-style-type: none"> Plotting Nyquist and Bode plots with statistical fitting Using equivalent circuit models to extract R_s, R_{ct}, C_{dl}, Warburg elements, etc. Evaluating charge-transfer kinetics and electrode–electrolyte interfacial behavior Clear interpretation of kinetic and diffusion contributions to OER performance <p>Tafel Plotting & Mechanistic Insights</p> <ul style="list-style-type: none"> Generating accurate Tafel plots (η vs $\log j$) with linear fitting Extracting Tafel slopes and exchange current density (j_0) Identifying dominant OER pathway (Adsorption/Deprotonation, M-O intermediate formation, etc.) Benchmarking catalyst efficiency against literature values <p>Stability Study Evaluation</p> <ul style="list-style-type: none"> Plotting CA/CP curves for long-term durability evaluation Analysing current retention, potential drift, and stability failure modes Comparing before/after LSV curves to assess surface activation or degradation Interpretation of catalyst robustness for real-world operation <p>Comprehensive Reporting</p> <ul style="list-style-type: none"> Providing high-quality, publication-ready plots for all techniques Presenting extracted parameters (η, Tafel slope, R_{ct}, C_{dl}, etc.) in structured tables Delivering in-depth interpretation of catalyst activity, kinetics, mechanism, and durability Summarising key findings in a concise, manuscript-friendly format | | |

| | | | |
|---|---|--------|----|
| | <p>Importance of These Analyses for OER Research & Manuscripts</p> <ul style="list-style-type: none"> • These characterisations together establish the four essential pillars of OER evaluation: Activity, kinetics, mechanism, and stability. • They provide the complete electrochemical evidence required for drafting a scientifically strong, publication-ready OER manuscript. • Reviewers expect LSV, CV, EIS, Tafel, and stability analyses for validating OER performance, benchmarking with existing catalysts, and demonstrating reliability and reproducibility. | | |
| 3 | Supercapacitor (CV, GCD, EIS, Stability)-3 electrode (Half Cell) | ₹14999 | ** |
| | <p>Analysis & Reporting Includes:**</p> <p>CV (Cyclic Voltammetry) Analysis</p> <ul style="list-style-type: none"> • Plotting CV curves at multiple scan rates for electrochemical profiling • Calculating specific capacitance (F/g, F/cm², or F/cm³) from CV curves • Identifying redox peaks (if pseudocapacitive) and differentiating EDLC vs. pseudocapacitive behaviour • Assessing rate capability through shape retention at higher scan rates <p>GCD (Galvanostatic Charge–Discharge) Analysis</p> <ul style="list-style-type: none"> • Plotting GCD curves at various current densities • Extracting specific capacitance, energy density, and power density • Determining IR-drop and interpreting internal resistance contribution • Evaluating charge–discharge symmetry to understand reversibility <p>EIS (Electrochemical Impedance Spectroscopy) Modelling</p> <ul style="list-style-type: none"> • Plotting Nyquist and Bode plots with accurate fitting • Extracting R_s (solution resistance), R_{ct} (charge-transfer resistance), and Warburg elements • Interpreting ion diffusion, interface behaviour, and electrode kinetics • Modelling capacitance response using appropriate equivalent circuits <p>Stability / Cycling Performance</p> <ul style="list-style-type: none"> • Plotting capacitance retention over extended cycles (CV/GCD) • Evaluating coulombic efficiency and charge–discharge stability • Comparing pre- and post-cycling performance to assess material durability • Identifying degradation mechanisms through curve distortion and IR-drop changes <p>Comprehensive Reporting</p> <ul style="list-style-type: none"> • Providing clean, publication-ready plots for CV, GCD, EIS, and cycling data • Presenting extracted parameters (specific capacitance, R_s/R_{ct}, retention %, IR-drop, etc.) in structured tables | | |

| | | | |
|---|--|-------|----|
| | <ul style="list-style-type: none"> Delivering clear interpretation of capacitive behaviour, kinetics, and long-term performance Summarising results in a concise, manuscript-friendly technical report <p>Importance of These Analyses for Supercapacitor Research & Manuscripts</p> <ul style="list-style-type: none"> Together, CV, GCD, EIS, and cycling stability provide the complete electrochemical fingerprint of a supercapacitor material. These analyses validate charge storage mechanism, performance, efficiency, kinetics, and durability, all of which are essential for a high-quality, publication-ready manuscript. Reviewers expect these four core techniques for benchmarking new electrode materials against state-of-the-art supercapacitors. | | |
| 4 | Supercapacitor (CV, GCD, EIS, Stability)-2 electrode (Full cell) | ₹9999 | ** |
| | <p>Analysis & Reporting Includes:**</p> <p>CV (Cyclic Voltammetry) Analysis</p> <ul style="list-style-type: none"> Plotting CV curves at multiple scan rates for full-cell performance evaluation Calculating specific capacitance (F/g, F/cm², or F/cm³) from CV curves Assessing device behaviour (EDLC vs. pseudocapacitive characteristics) Evaluating rate capability through CV shape retention at higher scan rates Determining voltage window suitability and operational stability <p>GCD (Galvanostatic Charge–Discharge) Analysis</p> <ul style="list-style-type: none"> Plotting charge–discharge curves at several current densities Extracting: <ul style="list-style-type: none"> Specific capacitance of the full cell Energy density (Wh/kg or Wh/L) Power density (W/kg or W/L) Measuring IR-drop to interpret internal resistance and electrode balancing Identifying charge–discharge symmetry for reversibility assessment <p>EIS (Electrochemical Impedance Spectroscopy) Modelling</p> <ul style="list-style-type: none"> Plotting Nyquist and Bode plots with high-quality curve fitting Extracting R_s (solution resistance), R_{ct} (charge-transfer resistance), and diffusion-related components Analysing ion transport efficiency and electrode–electrolyte interactions Determining device-level impedance and performance-limiting parameters Equivalent circuit modelling for full-cell devices <p>Stability / Cycling Performance</p> <ul style="list-style-type: none"> Plotting long-term cycling data (CV/GCD) for capacitance retention Evaluating coulombic efficiency, charge–discharge durability, and performance degradation | | |

- Comparing pre- and post-cycling CV/GCD curves to identify changes in behaviour
- Analysing loss mechanisms such as electrode degradation, electrolyte breakdown, and resistance increase

Comprehensive Reporting

- Providing clean, publication-ready graphs for CV, GCD, EIS, and cycling stability
- Presenting key parameters in tables:
 - Specific capacitance
 - Energy & power density
 - R_s/R_{ct} values
 - Retention %, coulombic efficiency
 - IR-drop trend analysis
- Delivering a concise, interpretation-rich report suitable for journal manuscripts, theses, or client presentations
- Highlighting comparative performance with literature benchmarks

Instructions

1. The prices listed are **indicative** and provided **for reference only**.
2. **Final charges** will be communicated after further **email or phone discussion**, based on a clear understanding of the **actual testing and consultancy requirements**.
3. Prices are **subject to change** depending on **sample nature, testing complexity, and scope of work**.
4. **Testing duration and slot allocation** depend on **facility and instrument availability**.
5. The testing timeline may vary from **1 to 3 weeks after receipt of payment**.
6. All the above prices are exclusive of GST. **An additional 18% GST will be applicable** and added to the final invoice as per Government of India regulations.

For ESSNANO MATERIALS LAB Pvt. Ltd.



Director